


C A G E N I X™

ACCUFRA[®]ME IC™

CUSTOM CERAMIC CROWNS ON OUR ACCUFRA[®]ME FRAMEWORK FOR ABOUT THE SAME PRICE AS A TRADITIONAL FIXED HYBRID.



HOW DOES ACCUFRA[®]ME IC STACK UP VS

FIXED HYBRIDS

SUPERIOR AESTHETICS

Ceramic crowns available in zirconia or e.Max with optional cutbacks for porcelain stacking

SUPERIOR DURABILITY

Custom ceramic crowns perfectly mated to our award-winning AccuFrame Plus Titanium frameworks

SUPERIOR FUNCTION

Excellent transferral of occlusal loads improves chewing performance as well as gingival durability

SUPERIOR REPAIRABILITY

Quick and Easy crown replacement via our lingual retrieval notch and digitally stored crown files

INDIVIDUAL CROWNS

SUPERIOR VALUE

AccuFrame IC pricing offers savings of up to 60% over traditional individual crown restorations

SUPERIOR PROCESSING

Easy system for processing acrylic gingiva means savings both in time and money

SUPERIOR TURN-AROUND

Fewer steps and patient visits plus quicker manufacturing turnaround means lower costs and quicker deliverability

SUPERIOR REPAIRABILITY

Our patent-pending crown and post design significantly improves ability to remove/repair the restoration.

For more information call us at (866) 964-5736 or visit cagenix.com

ACCUFAME[®] IC

INDIVIDUAL CROWN SYSTEM

Introducing AccuFrame IC

For years, the only options for fixed-hybrid cases have been traditional restorations with common denture teeth that are affordable but lack durability, or extremely expensive, complex individual crown options. Our AccuFrame IC (Individual Crown) System is removing the price barrier to improved look and performance of your fixed-hybrid cases — and making it possible for you to offer the highest quality full-arch restorations at truly affordable prices.

AccuFrame IC restorations utilize individual crowns in a choice of ceramic materials, including high-strength zirconia, e.max[®] or even layered porcelain, on our award-winning AccuFrame or AccuFrame Plus titanium frameworks. The result? Guaranteed fit, superior appearance and long-lasting performance.



AccuFrame IC vs Other Options

As you can see below, other restorations lack a means for the denture teeth to directly engage the bone structure. AccuFrame IC's titanium pontic form and framework structure has been specifically engineered to mate perfectly with the custom crowns, thereby transferring the occlusal force directly to the bone. This results in the best-performing, longest-lasting restoration available.

What type of patients benefit from Individual Crowns

- Patients seeking most realistic appearance available
- Patients seeking to minimize return visits for repair/rework
- Patients seeking to regain full functionality of all teeth
- Patients with increased wear traits (like Bruxers or grinders)



Traditional Denture
No structural contact between occlusal surface and bone



Fixed Hybrid
Frame engages bone but teeth supported via acrylic gingiva only



AccuFrame IC
Direct transfer of occlusal loads to bone for optimal performance

C A G E N I X[™]

For more info call (866) 964-5736 or visit www.cagenix.com

v0214.01